

## PROM2 Antibody

Catalog No: #62209

Package Size: #62209 100ul

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | PROM2 Antibody                                                                                                        |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB                                                                                                                    |
| Species Reactivity    | Human; Mouse; Rat                                                                                                     |
| Specificity           | The antibody detects endogenous level of total PROM2 protein                                                          |
| Immunogen Type        | Peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from part region of human protein                                                         |
| Target Name           | PROM2; PROML2; UNQ2521/PRO6014;                                                                                       |
| Calculated MW         | Predicted band size: 91 kDa                                                                                           |
| Formulation           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                        |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                       |

## Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.