

ECM29 Antibody

Catalog No: #62212

Package Size: #62212 100ul

Description

Product Name	ECM29 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human; Mouse
Specificity	The antibody detects endogenous level of total ECM29 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human ECM29 AA range: 1599-1649
Conjugates	Unconjugated
Target Name	ECM29
Calculated MW	203kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.