

C6orf173 Antibody

Catalog No: #62219

Package Size: #62219 100ul

Description

Product Name	C6orf173 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Applications	WB
Species Reactivity	Human, Mouse, Rat
Specificity	The antibody detects endogenous level of total C6orf173 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human C6orf173
Conjugates	Unconjugated.
Target Name	CENPW
Calculated MW	10kD
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.