

M4A6A Antibody

Catalog No: #62220

Package Size: #62220 100ul


Description

Product Name	M4A6A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total M4A6A protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human M4A6A AA range: 8-58
Conjugates	Unconjugated.
Target Name	M4A6A
Calculated MW	27kD
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Images



All lanes : M4A6A Antibody at 1/1k dilution
Lane 1 : HEK293 whole cell lysates lysates
Lysates/proteins at 20 ug per lane.
Secondary
All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 27 kDa
Observed band size: 27 kDa
Exposure time: 7 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.