

CTC1 Antibody

Catalog No: #62224

Package Size: #62224 100ul

Description

Product Name	CTC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Applications	WB
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous level of total CTC1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human CTC1
Conjugates	Unconjugated.
Target Name	CTC1
Calculated MW	134kD
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.