

IFM5 Antibody

Catalog No: #62225

Package Size: #62225 100ul

Description

Product Name	IFM5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen
Applications	WB
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous level of total IFM5 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human IFM5 AA range: 1-51
Conjugates	Unconjugated.
Target Name	IFM5
Calculated MW	15kD
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.