

F131C antibody

Catalog No: #62226

Package Size: #62226 100ul

Description

Product Name	F131C antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of human F131C.
Immunogen Type	Peptide
Immunogen Description	Peptide sequence around aa.130-180 derived from human F131C.
Target Name	F131C
Calculated MW	31kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.