

ZNF558 antibody

Catalog No: #62236

Package Size: #62236 100ul


Description

Product Name	ZNF558 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of ZN558.
Immunogen Description	Synthesized peptide derived from human ZN558 AA range: 230-280.
Target Name	ZNF558
Calculated MW	44kD
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from U2OS cells using ZNF558 antibody at 1:1000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.