

LSM2 Rabbit pAb

Catalog No: #62239

Package Size: #62239-1 50ul #62239-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LSM2 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of human LSM2
Immunogen Description	Synthesized peptide derived from human LSM2
Target Name	LSM2 C6orf28 G7B
Other Names	U6 snRNA-associated Sm-like protein LSm2 ; Protein G7b ; Small nuclear ribonuclear protein D homolog ; snRNP core Sm-like protein Sm-x5
Calculated MW	11kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-1:2000; IHC 1:100-1:300; IF 1:200-1:1000; ELISA 1:5000

Images



Western blot analysis of lysates from Jurkat cells, primary antibody was diluted at 1:1000, 4°over night



Immunohistochemical analysis of paraffin-embedded human tonsil.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.