

# BCKDHA (Phospho Ser292/Ser337) Antibody

Catalog No: #62242

Package Size: #62242-1 100ul


## Description

|                       |                                                                                                                                                                |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | BCKDHA (Phospho Ser292/Ser337) Antibody                                                                                                                        |
| Host Species          | Rabbit                                                                                                                                                         |
| Clonality             | Polyclonal                                                                                                                                                     |
| Purification          | The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns. |
| Applications          | IHC                                                                                                                                                            |
| Species Reactivity    | Human; Mouse; Rat                                                                                                                                              |
| Specificity           | BCKDHA (Phospho Ser292/Ser337) Antibody detects endogenous levels of BCKDHA only when phosphorylated at Ser292/337.                                            |
| Immunogen Description | A synthesized peptide derived from human BCKDHA around the phosphorylation site of Ser292/337.                                                                 |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                                                                |

## Application Details

IHC 1:50-1:200

## Images



BCKDHA (Phospho Ser292/Ser337) Antibody at 1/100 staining Mouse heart tissue by IHC-P.  
The sample was formaldehyde fixed and a heat mediated antigen retrieval step in citrate buffer was performed.  
The sample was then blocked and incubated with the primary antibody at 4°C overnight. An HRP conjugated anti-Rabbit antibody was used as the secondary antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.