

LPS/Lipopolysaccharide Mouse mAb

Catalog No: #62256

Package Size: #62256-1 100ul



Description

Product Name	LPS/Lipopolysaccharide Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4520
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA,IF
Species Reactivity	Escherichia coli
Formulation	0.01M PBS, pH 7.4.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

ELISA 1:5000-100000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.