

LPP1 Rabbit pAb

Catalog No: #62263

Package Size: #62263-1 50ul #62263-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LPP1 Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of LPP1 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human LPP1
Conjugates	Unconjugated
Calculated MW	31kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-2000; IHC 1:50-300; ELISA 1:2000-20000

Images



Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4°over night



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.