

GPR35 Rabbit pAb

Catalog No: #62264

Package Size: #62264-1 50ul #62264-2 100ul

Description

Product Name	GPR35 Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	GPR35 Polyclonal Antibody detects endogenous levels of GPR35 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human GPR35.
Conjugates	Unconjugated
Target Name	GPR35
Other Names	GPR35 ; G-protein coupled receptor 35 ; Kynurenic acid receptor ; KYNA receptor
SDS-PAGE MW	34kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-1:2000; IF 1:200-1:1000; ELISA 1:20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.