

Osteocalcin Rabbit mAb

Catalog No: #62265

Package Size: #62265-1 50ul #62265-1 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Osteocalcin Rabbit mAb
Host Species	rabbit
Clonality	monoclonal
Isotype	IgG,Kappa
Purification	Protein A
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	BGLAP
Calculated MW	11kDa
SDS-PAGE MW	22kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year


Application Details

IHC 1:2000-1:10000; WB 1:100-1:500; IF 1:200-1:1000; ELISA 1:5000-1:20000;


Images




Various whole cell lysates were separated by 4-8% SDS-PAGE, and the membrane was blotted with anti-Osteocalcin antibody. The HRP-conjugated Goat anti-Rabbit IgG (H + L) antibody was used to detect the antibody. Lane 1: A549 Lane 2: HepG2 Lane 3: Mouse lung
Predicted band size: 11kDa Observed band size: 22kDa




Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night



Mouse liver was stained with anti-Osteocalcin rabbit antibody



Human hepatocellular carcinoma was stained with anti-Osteocalcin rabbit antibody



Immunofluorescence analysis of A549.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.