

EpCAM Rabbit mAb

Catalog No: #62272

Package Size: #62272-1 50ul #62272-2 100ul


Description

Product Name	EpCAM Rabbit mAb
Host Species	rabbit
Clonality	monoclonal
Isotype	IgG,Kappa
Purification	Protein A
Applications	WB;IHC;IF;IP;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	EPCAM GA733-2 M1S2 M4S1 MIC18 TACSTD1 TROP1
Other Names	EPCAM ; GA733-2 ; M1S2 ; M4S1 ; MIC18 ; TACSTD1 ; TROP1 ; Epithelial cell adhesion molecule ; Ep-CAM ; Adenocarcinoma-associated antigen ; Cell surface glycoprotein Trop-1 ; Epithelial cell surface antigen ; Epithelial glycoprotein ; EGP ; Epithelial gly
Calculated MW	35kDa
SDS-PAGE MW	40kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year


Application Details

WB 1:1000-1:5000; IHC 1:200-1:1000;IF 1:200-1:1000; ELISA 1:5000-1:20000; IP 1:50-1:200,


Images




Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-EpCAM antibody. The HRP-conjugated Goat anti-Rabbit IgG(H + L) antibody was used to detect the antibody. Lane 1: HCT-116 Lane 2: MCF7 Predicted band size: 35kDa Observed band size: 40kDa




Rat colon was stained with anti-EpCAM rabbit antibody




Human liver was stained with anti-EpCAM rabbit antibody




Mouse colon was stained with anti-EpCAM rabbit antibody



Rat prostate was stained with anti-EpCAM Rabbit antibody



Rat prostate was stained with anti-EpCAM Rabbit antibody



Mouse colon was stained with anti-EpCAM rabbit antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.