

EMC2 Rabbit pAb

Catalog No: #62275

Package Size: #62275-1 50ul #62275-2 100ul


Description

Product Name	EMC2 Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of EMC2 at Human/Mouse/Rat
Immunogen Description	Synthesized peptide derived from human EMC2
Conjugates	Unconjugated
Other Names	EMC2;KIAA0103;TTC35
Calculated MW	33kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.