

NDUAC Rabbit pAb

Catalog No: #62277

Package Size: #62277-1 50ul #62277-2 100ul


Description

Product Name	NDUAC Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse
Specificity	NDUAC Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from human protein .
Conjugates	Unconjugated
Target Name	NDUFA12 DAP13
SDS-PAGE MW	15kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000; ELISA 1:5000-20000

Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.