

FoxO1A (Acetyl Lys245) Rabbit pAb

Catalog No: #62278

Package Size: #62278-1 50ul #62278-2 100ul


Description

Product Name	FoxO1A (Acetyl Lys245) Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat FoxO1A (Acetyl Lys245).
Immunogen Description	Synthesized peptide derived from human FoxO1A (Acetyl Lys245)
Conjugates	Unconjugated
Target Name	FOXO1 FKHR FOXO1A
Other Names	Forkhead box protein O1 ; Forkhead box protein O1A ; Forkhead in rhabdomyosarcoma
SDS-PAGE MW	72kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000-2000; ELISA 1:5000-20000

Images



Western Blot analysis of HEK 293 cells, cell treated or untreated by TSA 400nM 24h.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.