

CD28 Rabbit mAb

Catalog No: #62279

Package Size: #62279-1 50ul #62279-2 100ul


Description

Product Name	CD28 Rabbit mAb
Host Species	rabbit
Clonality	monoclonal
Isotype	IgG,Kappa
Purification	Protein A
Applications	WB;IHC;IF;IP;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	CD28
Other Names	CD28 ; T-cell-specific surface glycoprotein CD28 ; TP44 ; CD28
Calculated MW	25kDa
SDS-PAGE MW	50kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year


Application Details

WB 1:2000-1:10000; IHC 1:200-1:1000; IF 1:200-1:1000; IP 1:50-1:200; ELISA 1:5000-1:20000

Images



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-CD28 antibody.



Human tonsil was stained with anti-CD28 Rabbit antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.