

Cdc6 Rabbit pAb

Catalog No: #62280

Package Size: #62280-1 50ul #62280-2 100ul


Description

Product Name	Cdc6 Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	Cdc6 Polyclonal Antibody detects endogenous levels of Cdc6 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human CDC6.
Conjugates	Unconjugated
Target Name	CDC6
Other Names	CDC6 ; CDC18L ; Cell division control protein 6 homolog ; CDC6-related protein ; Cdc18-related protein ; HsCdc18 ; p62 ; cdc6 ; HsCDC6
SDS-PAGE MW	60kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-1:2000; IHC 1:100-1:300; ELISA 1:20000; IF 1:50-200

Images



Western blot analysis of lysates from K562 cells, using CDC6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using CDH9 antibody.



Immunohistochemical analysis of paraffin-embedded human small intestinal carcinoma tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.