

cPLA2 Rabbit mAb

Catalog No: #62289

Package Size: #62289-1 48T #62289-2 96T


Description

Product Name	cPLA2 Rabbit mAb
Host Species	rabbit
Clonality	monoclonal
Isotype	IgG,Kappa
Purification	Protein A
Applications	WB;IF;IP;ELISA
Species Reactivity	Human;Mouse
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	PLA2G4A
Other Names	PLA2G4A ; CPLA2 ; PLA2G4 ; Cytosolic phospholipase A2 ; cPLA2 ; Phospholipase A2 group IVA
Calculated MW	85kDa
SDS-PAGE MW	95kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

WB 1:2000-1:10000; IF 1:200-1:1000; IP 1:50-1:200; ELISA 1:5000-1:20000

Images



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-cPLA2 antibody. The HRP-conjugated Goat anti-Rabbit IgG (H + L) antibody was used to detect the antibody. Lane 1: HEK293 Lane 2: NIH-3T3 Predicted band size: 85kDa Observed band size: 95kDa

Note: This product is for in vitro research use only and is not intended for use in humans or animals.