

Dectin-2 Rabbit pAb

Catalog No: #62291

Package Size: #62291-1 50ul #62291-2 100ul


Description

Product Name	Dectin-2 Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Dectin-2 Polyclonal Antibody detects endogenous levels of Dectin-2 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Dectin-2.
Conjugates	Unconjugated
Target Name	CLEC6A
Other Names	CLEC6A ; CLECSF10 ; DECTIN2 ; C-type lectin domain family 6 member A ; C-type lectin superfamily member 10 ; Dendritic cell-associated C-type lectin 2 ; DC-associated C-type lectin 2 ; Dectin-2
SDS-PAGE MW	23kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500-1:2000; IHC: 1:100-1:300; ELISA 1:10000; IF 1:50-200


Images



Western blot analysis of SH-SY5Y 293T using Dectin-2 antibody. Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:20000



Western blot analysis of K562 using Dectin-2 antibody. Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-heart, antibody was diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.