

GPER1 Rabbit mAb

Catalog No: #62292

Package Size: #62292-1 50ul #62292-2 100ul

Description

Product Name	GPER1 Rabbit mAb
Host Species	rabbit
Clonality	monoclonal
Isotype	IgG,Kappa
Purification	Protein A
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Conjugates	Unconjugated
Target Name	GPER CEPR CMKRL2 DRY12 GPR30
Other Names	G-protein coupled estrogen receptor 1 ; Chemoattractant receptor-like 2 ; Flow-induced endothelial G-protein coupled receptor 1 ; FEG-1 ; G-protein coupled receptor 30 ; GPCR-BR ; IL8-related receptor DRY12 ; Lymphocyte-derived G-protein coupled receptor ; LYGPR ; Membrane estrogen receptor ; mER ;
Calculated MW	42kDa
SDS-PAGE MW	55kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

IHC 1:1000-1:5000; WB 1:2000-1:10000; IF 1:200-1:1000; ELISA 1:5000-1:20000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.