

IGS11 Rabbit pAb

Catalog No: #62293

Package Size: #62293-1 50ul #62293-2 100ul


Description

Product Name	IGS11 Rabbit pAb
Host Species	rabbit
Clonality	polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	IGS11 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from human protein
Conjugates	Unconjugated
Target Name	IGSF11
Other Names	IGSF11;BTIGSF;CXADRL1;VSIG3
SDS-PAGE MW	47kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000; ELISA 1:5000-20000

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.