

Bcl-6 Mouse mAb

Catalog No: #62295

Package Size: #62295-1 50ul #62295-2 100ul


Description

Product Name	Bcl-6 Mouse mAb
Host Species	mouse
Clonality	monoclonal
Isotype	IgG3,Kappa
Purification	Protein G
Applications	WB;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of Bcl-6.
Immunogen Description	Synthesized peptide derived from human BCL6.
Conjugates	Unconjugated
Target Name	BCL6
Other Names	BCL6 ; BCL5 ; LAZ3 ; ZBTB27 ; ZNF51 ; B-cell lymphoma 6 protein ; BCL-6 ; B-cell lymphoma 5 protein ; BCL-5 ; Protein LAZ-3 ; Zinc finger and BTB domain-containing protein 27 ; Zinc finger protein 51
Calculated MW	78kDa
SDS-PAGE MW	78kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000; IF 1:100-500; ELISA 1:1000-5000

Images



Various whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Bcl-6 (PTR2017) antibody. The HRP-conjugated Goat anti-Mouse IgG(H + L) antibody was used to detect the antibody. Lane 1: Daudi Lane 2: Raji Lane 3: Ramos

Note: This product is for in vitro research use only and is not intended for use in humans or animals.