

ABCC4 Mouse mAb

Catalog No: #63004

Package Size: #63004-1 50ul #63004-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ABCC4 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human ABCC4 (AA: 631-692) expressed in E. Coli.
Target Name	ABCC4
Other Names	MRP4; MOATB; MOAT-B; EST170205
Accession No.	O15439
Calculated MW	150kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Flow cytometric analysis of A549 cells using ABCC4 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human endometrial cancer tissues using ABCC4 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human bladder cancer tissues using ABCC4 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.