

ACHE Mouse mAb

Catalog No: #63010

Package Size: #63010-1 50ul #63010-2 100ul


Description

Product Name	ACHE Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human,Mouse,Monkey,Rat
Immunogen Description	Purified recombinant fragment of human ACHE expressed in E. Coli.
Target Name	ACHE
Other Names	YT; ACEE; ARACHE; N-ACHE
Accession No.	P22303
Calculated MW	67.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400


Images




Western blot analysis using ACHE mouse mAb against PC-12 (1), Hela (2), mouse brain (3), NIH/3T3 (4), COS7 (5), Jurkat (6) and Raji (7) cell lysate.



Flow cytometric analysis of NIH/3T3 cells using ACHE mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissues using ACHE mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using ACHE mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.