

CD49E Mouse mAb

Catalog No: #63401

Package Size: #63401-1 50ul #63401-2 100ul


Description

Product Name	CD49E Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD49E (AA: extra 111-253) expressed in E. Coli.
Target Name	CD49E
Other Names	ITGA5; FNRA; VLA-5; VLA5A
Accession No.	P08648
Calculated MW	114.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD49E mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.