

CD57 Mouse mAb

Catalog No: #63406

Package Size: #63406-1 50ul #63406-2 100ul


Description

Product Name	CD57 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;FC
Species Reactivity	Human,Mouse,Rat,Monkey
Immunogen Description	Purified recombinant fragment of human CD57 (AA: 35-191) expressed in E. Coli.
Target Name	CD57
Other Names	B3GAT1; NK1; HNK1; LEU7; NK-1; GLCATP; GLCUATP
Accession No.	Q9P2W7
Calculated MW	38.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CD57 mouse mAb against SK-N-SH (1), SH-SY5Y (2), C6 (3), HL-60 (4), COS7 (5), Hela (6), and NIH/3T3 (7) cell lysate.



Flow cytometric analysis of HL-60 cells using CD57 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.