

CD57 Mouse mAb

Catalog No: #63407

Package Size: #63407-1 50ul #63407-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD57 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD57 (AA: 28-334) expressed in E. Coli.
Target Name	CD57
Other Names	NK1; CD57; HNK1; LEU7; NK-1; GLCATP; GLCUATP
Accession No.	Q9P2W7
Calculated MW	38.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200-1/1000FC:1/200-1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of Jurkat cells using CD57 mouse mAb (green) and negative control (red).



Flow cytometric analysis of HeLa cells using CD57 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human brain tissues using CD57 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.