

CD62P Mouse mAb

Catalog No: #63410

Package Size: #63410-1 50ul #63410-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CD62P Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD62P (AA: extra 42-208) expressed in E. Coli.
Target Name	CD62P
Other Names	SELP; CD62; GRMP; PSEL; GMP140; LECAM3; PADGEM
Accession No.	P16109
Calculated MW	90.8kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


FC:1/200 - 1/400ELISA:1/10000

Images

O.D. ELISA Result



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of THP-1 cells using CD62P mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.