

CD66A Mouse mAb

Catalog No: #63413

Package Size: #63413-1 50ul #63413-2 100ul


Description

Product Name	CD66A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD66A (AA: extra 65-201) expressed in E. Coli.
Target Name	CD66A
Other Names	CEACAM1; BGP; BGP1; BGPI
Accession No.	P13688
Calculated MW	57.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using CD66A mouse mAb against HepG2 (1), SK-MES-1 (2), Hela (3), and A549 (4) cell lysate.



Flow cytometric analysis of Hela cells using CD66A mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissues using CD66A mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD66A mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.