

CD69 Mouse mAb

Catalog No: #63414

Package Size: #63414-1 50ul #63414-2 100ul


Description

Product Name	CD69 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD69 expressed in E. Coli.
Target Name	CD69
Other Names	CLEC2C; CD69
Accession No.	Q07108
Calculated MW	22.5kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images



Western blot analysis using CD69 mouse mAb against, Jurkat (1), L1210 (2) and TPH-1 (3) cell lysate.



Flow cytometric analysis of Jurkat cells using CD69 mouse mAb (green) and negative control (purple).



Immunohistochemical analysis of paraffin-embedded human Tonsil tissues using anti-CD69 mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.