

CD79A Mouse mAb

Catalog No: #63417

Package Size: #63417-1 50ul #63417-2 100ul


Description

Product Name	CD79A Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD79A (AA: extra(33-143)) expressed in E. Coli.
Target Name	CD79A
Other Names	IGA; MB-1
Accession No.	P11912
Calculated MW	25kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using *** mouse mAb against Raji (1), Ramos (2), k562 (3), and MOLT4 (4) cell lysate.



Immunofluorescence analysis of HeLa cells using CD79A mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of Raji cells using CD79A mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.