

CD80 Mouse mAb

Catalog No: #63420

Package Size: #63420-1 50ul #63420-2 100ul


Description

Product Name	CD80 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of CD80 expressed in E. Coli.
Target Name	CD80
Other Names	CD28LG; LAB7; B7.1
Accession No.	P33681
Calculated MW	55kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000ELISA:1/10000

Images



Confocal Immunofluorescence analysis of BCBL-1 cells using anti-CD80 monoclonal antibody (green), showing membrane localization. Blue: DRAQ5 fluorescent DNA dye.



Immunohistochemical analysis of paraffin-embedded human Tonsil tissues using anti-CD80 mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.