

CD80 Mouse mAb

Catalog No: #63421

Package Size: #63421-1 50ul #63421-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD80 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human,Mouse,Rat,Rabbit
Immunogen Description	Purified recombinant fragment of human CD80 (AA: 35-288) expressed in E. Coli.
Target Name	CD80
Other Names	B7; BB1; B7-1; B7.1; LAB7; CD28LG; CD28LG1
Accession No.	P33681
Calculated MW	33kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/100 - 1/500FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of HL-60 cells using CD80 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD80 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rat spleen using CD80 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded Rabbit spleen using CD80 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.