

CD82 Mouse mAb

Catalog No: #63423

Package Size: #63423-1 50ul #63423-2 100ul


Description

Product Name	CD82 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	IHC;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human CD82 (AA: extra 111-228) expressed in E. Coli.
Target Name	CD82
Other Names	R2; 4F9; C33; IA4; ST6; GR15; KAI1; SAR2; TSPAN27
Accession No.	P27701
Calculated MW	29.6kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Flow cytometric analysis of NIH/3T3 cells using CD82 mouse mAb (green) and negative control (red).



Flow cytometric analysis of Hela cells using CD82 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD82 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human lymph tissues using CD82 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.