

CD85G Mouse mAb

Catalog No: #63426

Package Size: #63426-1 50ul #63426-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD85G Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD85G (AA: extra 24-178)) expressed in E. Coli.
Target Name	CD85G
Other Names	LILRA4;ILT7
Accession No.	P59901
Calculated MW	55.2kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200-1/1000 FC: 1/200-1/400 ELISA: 1/10000


Images



Black line: Control Antigen (100 ng); Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng)



Flow cytometric analysis of THP-1 cells using CD85G mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human tonsil tissues using CD85G mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.