

CD86 Mouse mAb

Catalog No: #63429

Package Size: #63429-1 50ul #63429-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CD86 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD86 expressed in E. Coli.
Target Name	CD86
Other Names	B70; B7-2; LAB72; CD28LG2; MGC34413; CD86
Accession No.	P42081
Calculated MW	52kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Western blot analysis using CD86 mouse mAb against L1210 (1) and MOLT-4 (2) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.