

CD88 Mouse mAb

Catalog No: #63431

Package Size: #63431-1 50ul #63431-2 100ul


Description

Product Name	CD88 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Applications	WB;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human CD88 (AA: extra mix) expressed in E. Coli.
Target Name	CD88
Other Names	C5AR1;C5A; C5AR; C5R1
Accession No.	P21730
Calculated MW	39kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000FC:1/200 - 1/400ELISA:1/10000


Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CD88 mouse mAb against mouse serum (1) tissue lysate.



Flow cytometric analysis of THP-1 cells using CD88 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.