

CD89 Mouse mAb

Catalog No: #63432

Package Size: #63432-1 50ul #63432-2 100ul


Description

Product Name	CD89 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD89 (AA: extra 22-227) expressed in E. Coli.
Target Name	CD89
Other Names	FCAR; FcalphaRI; CTB-61M7.2
Accession No.	P24071
Calculated MW	32.3kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of HL-60 cells using CD89 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.