

CD93 Mouse mAb

Catalog No: #63434

Package Size: #63434-1 50ul #63434-2 100ul


Description

Product Name	CD93 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IHC;IF;FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CD93 (AA: 474-535) expressed in E. Coli.
Target Name	CD93
Other Names	C1QR1; C1qRP; CDw93; ECSM3; MXRA4; C1qR(P); dJ737E23.1
Accession No.	Q9NPY3
Calculated MW	68.6kDa
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunofluorescence analysis of HeLa cells using CD93 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.



Flow cytometric analysis of HepG2 cells using CD93 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CD93 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.