

CDC25C Mouse mAb

Catalog No: #63440

Package Size: #63440-1 50ul #63440-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CDC25C Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB;IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CDC25C expressed in E. Coli.
Target Name	CDC25C
Other Names	CDC25; CDC25C
Accession No.	P30307
Calculated MW	58kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ELISA:1/10000

Images



Western blot analysis using anti-CDC25C mAb against Hela (1), K562 (2), PC-3 (3), HEK293 (4) and Raw264.7 (5) cell lysate.



Immunohistochemical analysis of paraffin-embedded human Testis tissues using anti-CDC25C mouse mAb

Note: This product is for in vitro research use only and is not intended for use in humans or animals.