

## CDH17 Mouse mAb

Catalog No: #63444

Package Size: #63444-1 50ul #63444-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                          |
|-----------------------|------------------------------------------------------------------------------------------|
| Product Name          | CDH17 Mouse mAb                                                                          |
| Host Species          | Mouse                                                                                    |
| Clonality             | Monoclonal                                                                               |
| Isotype               | Mouse IgG1                                                                               |
| Applications          | IHC                                                                                      |
| Species Reactivity    | Human                                                                                    |
| Immunogen Description | Purified recombinant fragment of human CDH17 (AA: extra(600-707)) expressed in E. Coli.  |
| Target Name           | CDH17                                                                                    |
| Other Names           | HPT1; CDH16; HPT-1                                                                       |
| Accession No.         | Q12864                                                                                   |
| Calculated MW         | 92.2kDa                                                                                  |
| Formulation           | Purified antibody in PBS with 0.05% sodium azide.                                        |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |


## Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000


## Images



Black line: Control Antigen (100 ng); Purple line: Antigen(10ng); Blue line: Antigen (50 ng); Red line: Antigen (100 ng);



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues using CDH17 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human colon tissues using CDH17 mouse mAb with DAB staining.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.