

CEA Mouse mAb

Catalog No: #63454

Package Size: #63454-1 50ul #63454-2 100ul


Description

Product Name	CEA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	IF
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CEA expressed in E. Coli.
Target Name	CEA
Other Names	CEA; CD66e; DKFZp781M2392; CEACAM5
Accession No.	P06731
Calculated MW	77kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

ICC:1/200 - 1/1000ELISA:1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Immunofluorescence analysis of PANC-1 cells using CEA mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.