

CEA Mouse mAb

Catalog No: #63455

Package Size: #63455-1 50ul #63455-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CEA Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4548
Isotype	Mouse IgG1
Applications	ELISA, IHC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CEA expressed in E. Coli.
Target Name	CEA
Other Names	CEA; CD66e; DKFZp781M2392; CEACAM5
Accession No.	P06731
Calculated MW	77kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IHC: 1/200 - 1/1000 ELISA: 1/10000

Images



Red: Control Antigen (100ng); Purple: Antigen (10ng); Green: Antigen (50ng); Blue: Antigen (100ng);



Immunohistochemical analysis of paraffin-embedded human rectum cancer tissues (left) and stomach cancer tissues (right) using CEA mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.