

CFHR5 Mouse mAb

Catalog No: #63462

Package Size: #63462-1 50ul #63462-2 100ul


Description

Product Name	CFHR5 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	WB;IHC;IF;FC
Species Reactivity	Human,Mouse
Immunogen Description	Purified recombinant fragment of human CFHR5 (AA: 344-569) expressed in E. Coli.
Target Name	CFHR5
Other Names	FHR5; CFHL5; FHR-5; CFHR5D
Accession No.	Q9BXR6
Calculated MW	64.4kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000IHC:1/200 - 1/1000ICC:1/200 - 1/1000FC:1/200 - 1/400ELISA:1/10000


Images




Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)




Western blot analysis using CFHR5 mouse mAb against HepG2 (1), K562 (2), L-02 (3), SK-Hep-1 (4), SMMC-7721 (5), and NIH/3T3 (6) cell lysate.




Immunofluorescence analysis of HeLa cells using CFHR5 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor- 555 phalloidin. Secondary antibody from Fisher (Cat#: 35503)



Flow cytometric analysis of HeLa cells using CFHR5 mouse mAb (green) and negative control (red).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissues using CFHR5 mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded human ovarian cancer tissues using CFHR5 mouse mAb with DAB staining.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.