

CHK1 Mouse mAb

Catalog No: #63464

Package Size: #63464-1 50ul #63464-2 100ul


Description

Product Name	CHK1 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Purified recombinant fragment of CHK1 expressed in E. Coli.
Target Name	CHK1
Other Names	CHEK1
Accession No.	O14757
Calculated MW	56kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500 - 1/2000 ELISA: 1/10000

Images



Western blot analysis using CHK1 mouse mAb against A431 (1), HeLa (2), NIH/3T3 (3) and K562 (4) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.