

CHRNA2 Mouse mAb

Catalog No: #63469

Package Size: #63469-1 50ul #63469-2 100ul


Description

Product Name	CHRNA2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Purified recombinant fragment of human CHRNA2 (AA: extra 27-264) expressed in E. Coli.
Target Name	CHRNA2
Accession No.	Q15822
Calculated MW	60kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB:1/500 - 1/2000ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Western blot analysis using CHRNA2 mouse mAb against SK-N-SH (1), SH-SY5Y (2), membrane protein of C6 (3), and SW480 (4) cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.