

CHRNA2 Mouse mAb

Catalog No: #63474

Package Size: #63474-1 50ul #63474-2 100ul


Description

Product Name	CHRNA2 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Applications	FC
Species Reactivity	Human
Immunogen Description	Purified recombinant fragment of human CHRNA2 (AA: extra 26-233) expressed in E. Coli.
Target Name	CHRNA2
Other Names	EFNL3; nAChRB2
Accession No.	P17787
Calculated MW	57kDa
Formulation	Purified antibody in PBS with 0.05% sodium azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

FC:1/200 - 1/400ELISA:1/10000

Images



Black line: Control Antigen (100 ng);Purple line: Antigen (10ng); Blue line: Antigen (50 ng); Red line:Antigen (100 ng)



Flow cytometric analysis of SH-SY5Y cells using CHRNA2 mouse mAb (green) and negative control (red).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.